



## TECHNICAL DETAILS



# Digital GSM Balance i9™ (ECO)

Economical Digital GSM Balance

Page 1 of 2

### UNMATCHED LIFETIME BUYBACK WARRANTY



### OVERVIEW



A Digital GSM Balance is used to determine GSM of all type of textile fabric, paper, plastic film, Leather, Paper & Packaging Industry etc. also to weigh different laboratory samples up to 600 gm with accuracy of 0.01 gm, so it is very useful to weigh samples & chemicals in textile laboratories, College & Universities.

**Paramount Digital GSM BALANCE i9™ (ECO MODEL)** is an advanced but Economical GSM calculating balance with the accuracy up to two decimal numbers. Users friendly Bigger Display and pan size. A.C adaptor is supplied with the equipment to run on 220V 50Hz A.C supply. A battery backup is also provided with it.

### SALIENT FEATURES



- External calibration
- Auto zero tracking
- Low battery indication
- Large weighing pan
- Net weight stable indication
- Large LCD display
- Kg, gm, oz, lb, tl, ct conversion
- Net weight /Stability indication

 **Paramount**<sup>®</sup>  
S I N C E 1 9 6 4  
YOUR GROWTH. OUR PASSION™

**Sales (24x7):** +91 999 999 1080, +91 9999 8899 33  
**Customer Care:** +91 999 999 8037, +91 999 00 999 61

[sales@paramountinstruments.com](mailto:sales@paramountinstruments.com)  
[www.paramountinstruments.com](http://www.paramountinstruments.com)



### TECHNICAL SPECIFICATIONS



Capacity of the Balance	: 600 gm.
Least Count	: 0.1 mg.
Platter size	: 15mm
Power source	: 1.5v x 2 xaa size batteries.
Ac Adapter	: 3v/100mA
Net Weight of the unit	: 617GM
Operating Temp.	: 0 to 40°C

### RELATED STANDARDS



Made As Per Customer's Requirements

### STANDARD ACCESSORIES



1. Main Unit.	: 1 No.
2. Acrylic Cover	: 1 No.
3. A.C. Adaptor	: 1 No.
4. Pan cover	: 1 No.
5. Calibration Certificate (Traceable to NPL)	: 1 No.
6. Inspection and Conformance Certificate	: 1 No.

### EXCLUSIVE BENEFITS



**WOW** Customer SERVICE



Complete **Peace of Mind**



**LIFETIME** Buy Back  
Warranty



ISO 9001 : 2015

QUALITY EXCELLENCE  
is Our Forte.



ISO 14001 : 2015

We CARE for Our  
Environment.



ISO 45001 : 2015

HEALTH & SAFETY  
Is Our Priority.



**Sales (24x7) :** +91 999 999 1080, +91 9999 8899 33  
**Customer Care :** +91 999 999 8037, +91 999 00 999 61

[sales@paramountinstruments.com](mailto:sales@paramountinstruments.com)  
[www.paramountinstruments.com](http://www.paramountinstruments.com)