



TECHNICAL SPECIFICATIONS



Size of the Sample to be cut : 152 x 13mm
 Overall Dimensions of the unit : 125mm x 530mm x 430mm
 5" x 21" x 17"
 Net Weight of the unit : 22 Kg. (48 lbs)

RELATED STANDARDS



TAPPI T 811 | IS 7063 | FEFCO TM 8

STANDARD ACCESSORIES



- 1. Main Unit : 01 No.
- 2. The User's Manual : 01 No.
- 3. Calibration Certificate : 01 No
- 4. Inspection & Conformance Certificate : 01 No

EXCLUSIVE BENEFITS



WOW Customer SERVICE



Complete **Peace of Mind**



LIFETIME Buy Back
 Warranty



QUALITY EXCELLENCE
 is Our Forte.

ISO 9001 : 2015



We CARE for Our
 Environment.

ISO 14001 : 2015



HEALTH & SAFETY
 Is Our Priority.

ISO 45001 : 2015



Sales (24x7) : +91 999 999 1080, +91 9999 8899 33
Customer Care : +91 999 999 8037, +91 999 00 999 61

sales@paramountinstruments.com

www.paramountinstruments.com